

Perfect Liquid Handling

PIPETTING WITH EASE AND COMFORT

NEW!



VITLAB 
Competence in Labware

Pipet controller for graduated, volumetric and blow-out pipettes

VITLAB maneus pipet controller

VITLAB specialises in accurate dosing, titration and pipetting and offers a wide range of highly-developed precision instruments for various Liquid Handling applications.

The new VITLAB maneus pipet controller allows both left and right-handed people to work easily, tirelessly

and safely with all common pipettes from 0.1 to 100 ml. Its easy handling combined with high sensitivity also lets inexperienced users set the meniscus precisely. The VITLAB maneus pipet controller is easy to take apart, very simple to clean and completely autoclavable.

