



Metanol, J.T. Baker®

Proveedor: Avantor

Sinónimos: Destilado metílico , Alcohol metílico

Formula: CH₃OH

MW: 32,04 g/mol

Número MDL:

MFCD00004595

Núm. CAS: 67-56-1

UN: 1230

ADR: 3,II



[Order Now](#)

Resultados test especificaciones

Exceeds ACS Specifications	
Assay (by GC) (corrected for water)	min. 99.8 %
Appearance	clear
Carbonyl Compounds (max. 0.001% each of Acetone, Formaldehyde and Acetaldehyde)	passes test
Color (APHA)	max. 10
Residue after Evaporation	max. 0.001 %
Solubility in Water	passes test
Substances Darkened by H ₂ SO ₄	passes test
Substances Reducing KMnO ₄	passes test
Titration Acid (meq/g)	max. 0.0003

Titration Base (meq/g)	max. 0.0002
Water (H ₂ O)	max. 0.01 %
Trace Impurities (in ppm)	
Aluminium (Al)	max. 0.5
Boron (B)	max. 0.02
Barium (Ba)	max. 0.1
Calcium (Ca)	max. 0.5
Cadmium (Cd)	max. 0.05
Cobalt (Co)	max. 0.02
Chromium (Cr)	max. 0.02
Copper (Cu)	max. 0.02
Iron (Fe)	max. 0.1
Magnesium (Mg)	max. 0.1
Manganese (Mn)	max. 0.02
Nickel (Ni)	max. 0.02
Lead (Pb)	max. 0.1
Tin (Sn)	max. 0.1
Zinc (Zn)	max. 0.1

[PEDIDO](#)[DOCUMENTACIÓN](#)

Recommended for You